

Certificate of Decontamination

Instrument manufacturer & model: _____

Instrument serial number: _____

Has this equipment been exposed to any infectious agents? YES / NO

If YES, please indicate which agent and bio-safety hazard class:

Has this equipment been exposed to radioactive substances? YES / NO

If YES, please indicate isotopes and general quantities used:

Has this equipment been exposed to toxic or carcinogenic materials? YES / NO

If YES, please indicate which types:

Describe how the equipment was cleaned or disinfected by checking one or more of the decontamination procedures outlined below:

- Equipment not contaminated:** User acknowledges that instrumentation has not come into contact with, or contains, any infectious agents, radioactive materials, or toxic chemicals.
- Exposure to infectious agents:** The equipment was thoroughly decontaminated by spraying it or wiping it down with a freshly prepared solution of 10% bleach in water. The disinfectant remained on the equipment for a minimum of 15 minutes before rinsing with water.
- Exposure to radioactive materials:** The equipment was thoroughly decontaminated by cleaning it with a commercially available decontaminant (i.e. Radiacwash®) or equivalent. The instrument was subsequently surveyed via wipe testing using the appropriate instrumentation and the calibration report and results are attached to this form.
- Exposure to toxic chemicals:** The equipment was thoroughly decontaminated by rinsing areas exposed to chemicals with a solvent such as methanol or ethanol, then rinsing with detergent and water.

By signing this document, the following party warrants that all appropriate decontamination procedures as outlined above have been performed and that the equipment is completely decontaminated.

Company: _____

Signature: _____

Print Name: _____

Date: _____